



Figure A.2: Schematic diagrams of the high voltage controller showing the sections associated with controlling VF2, VF3, and VF4. The dashed rectangles enclose parts located on the main circuit board. The numbers and letters contained in boxes refer to pins on the circuit board edge connector, and the characters in quotes refer to labels on the front panel of the chassis.






(to bottom el. of X-Y 2) (to bottom el. of X-Y 1) (to right el. of X-Y 1) (to left el. of X-Y 1) (to right el. of X-Y 2) (to left el. of X-Y 2) (to top el. of X-Y 2) (to top el. of X-Y 1) X-Y 2 Right/Left switch X-Y 1 Right/Left switch X-Y 2 Right/Left switch X-Y 2 Up/Down switch X-Y 2 Up/Down switch X-Y 1 Up/Down switch X-Y 1 Right/Left switch Up/Down switch X-Y 1 C1-H; C1-F: C1-G C1-D: C1-E: C1-B: CIPC


Parts List


U1–U4: Burr-Brown 3582J U5–U8: LM343 U9: REF-01HP Unless otherwise noted, all resistors are 1/4 W, 1% metal film.


Notes

- 1. All three power supplies are separately fused, apart from the fuse in the line filter. These are labeled F1, F2, and F3 on the back panel. F1 goes to the HHV-150, F2 goes to the HHV-50, and F3 goes to the Power-One.
- 2. U1-U4 were originally 3584 op-amps, which require external compensation. Consequently, there are components attached to pin 8 of each 3582, which is not connected. These components can be ignored.









Parts List

U1, U2: Burr-Brown 3582J U3: OP-07DP U4: REF-01HP Unless otherwise noted, all resistors are 1/4 W, 1% metal

Notes

1. All three power supplies are separately fused, apart from the fuse in the line filter. These fuses are labeled F1, F2, and F3 on the back panel. F1 and F2 go to the Ferrotran supplies, and F3 goes to the Power-One.




Figure A.8: The schematic diagram of the monochromator low voltage supply. The dashed rectangle encloses those components located on the circuit board. The numbers enclosed by squares refer to pins on the circuit board edge connector, and the characters in quotes are labels on the front or back panel.

Parts List

U1: OP-07DP
U2, U3: Linear Technology
LTH012CN
U4: REF-01HP
Unless otherwise noted, all
resistors are 1/4 W, 1% metal

Notes

The labels 'a', 'b', and 'c' on the circuit board denote test points.

a=output(+) + 400

b=output(+) + 4000

c= the inverting input of U2-nominally 0.V.

OV.

C= the inverting input of U2-nominally 0.V.

C= the inverting input of U2-nominally 0.V.

OV.

The labels 'a', 'b', and 'c' on the circuit board denote test points.

C= the inverting input of U2-nominally 0.V.

OV.

The labels 'a', 'b', and 'c' on the circuit board denote test points.

C= the inverting input of U2-nominally 0.V.

OV.

The labels 'a', 'b', and 'c' on the circuit board denote test points.

HAA15-0.8-A Power-One PMT-50A/N 0 to -5 kV ±15 V Bertan PMT-50A/P 0 to +5 kV Circuit Board Bertan 2007 2 2 3 **5** 5 阊 Corcom line filter monitor coutput (1V/1 keV or 1V/10 keV) program [voltage output

